

Catalog #: 52421

Aliases

CKAP4; Cytoskeleton Associated Protein 4; P63; Cytoskeleton-Associated Protein 4; CLIMP-63; ERGIC-63; CLIMP63; 63-KDa Cytoskeleton-Linking Membrane Protein; Transmembrane Protein (63kD), Endoplasmic Reticulum/Golgi Intermediate Compartment; Type-II Transmembrane Protein P63; 63 KDa Membrane Protein; Climp-63

Background

Gene Name: CKAP4

NCBI Gene Entry: [10970](#)

UniProt Entry: [Q07065](#)

Application Information

Molecular Weight: Predicted, 66 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human CKAP4

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

