

Anti-M13 bacteriophage Mouse Monoclonal Antibody



Catalog #: 52450

Background

M13 Phage; M13 Filamentous Bacteriophage; M13 Virus; Filamentous Coliphage M13

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 25GB50

Species Reactivity: Bacteriophage M13

Applications Tested: Western blotting (WB)

Immunogen

M13 bacteriophage

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS containing 0.09% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:200-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



7 μ L of M13 bacteriophage lysate was run on 6-18% SDS-PAGE under reducing conditions and blotted onto nitrocellulose membrane. M13 bacteriophage antibody (Cat#52450) at 1 μ g/mL was used as the primary antibody and peroxidase conjugated goat anti-mouse IgG was used as the secondary antibody. Image was developed using FeQ™ ECL Substrate Kit.