

Anti-Adenylate Kinase 9 Rabbit Polyclonal Antibody



Catalog #: 52454

Aliases

AK9; Adenylate Kinase 9; DJ70A9.1; C6orf199; C6orf224; AKD1; AKD2; Adenylate Kinase Domain Containing 2; Adenylate Kinase Domain Containing 1; FLJ42177; FLJ25791; MGC26954; Adenylate Kinase Domain-Containing Protein 1; Adenylate Kinase Domain-Containing Protein 2; Chromosome 6 Open Reading Frame 224; Chromosome 6 Open Reading Frame 199; EC 2.7.4.4; EC 2.7.4.6; SPGF89; AK 9

Background

Gene Name: AK9

NCBI Gene Entry: [221264](#)

UniProt Entry: [Q5TCS8](#)

Application Information

Molecular Weight: Predicted, 221 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human AKD2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

