

Catalog #: 52459

Aliases

NMNAT3; Nicotinamide Nucleotide Adenylyltransferase 3; PNAT3; Nicotinamide/Nicotinic Acid Mononucleotide Adenylyltransferase 3; Nicotinate-Nucleotide Adenylyltransferase 3; Pyridine Nucleotide Adenylyltransferase 3; NMN/NaMN Adenylyltransferase 3; NaMN Adenylyltransferase 3; NMN Adenylyltransferase 3; PNAT-3; Nicotinamide Mononucleotide Adenylyltransferase 3; Nicotinamide-Nucleotide Adenylyltransferase 3; EC 2.7.7.18; EC 2.7.7.1; FKSG76

Background

Gene Name: NMNAT3

NCBI Gene Entry: [349565](#)

UniProt Entry: [Q96T66](#)

Application Information

Molecular Weight: Predicted, 28 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human NMNA3

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

SUPPORT


SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Validation Data



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-NMNAT3 antibody (Cat#52459). Total lysates (30 µg) were loaded and separated by SDS-PAGE. The blot was incubated with anti-NMNAT3 antibody (Cat#52459, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).