

Catalog #: 5246

Aliases

RALGAPA1; Ral GTPase Activating Protein Catalytic Subunit Alpha 1; GRIPE; Tuberin-Like Protein 1; RalGAPalpha1; KIAA0884; Tulip1; GARNL1; Ral GTPase-Activating Protein Subunit Alpha-1; GTPase Activating Rap/RanGAP Domain-Like 1; GAP-Related Interacting Protein To E12; DKFZp667F074; P240; Ral GTPase Activating Protein, Alpha Subunit 1 (Catalytic); Ral GTPase Activating Protein Catalytic Alpha Subunit 1; GTPase-Activating Rap/Ran-GAP Domain-Like 1; GTPase Activating RANGAP Domain-Like 1; GAP-Related-Interacting Partner To E12; NEDHRIT; TULIP1

Background

Gene Name: RALGAPA1

NCBI Gene Entry: [253959](#)

UniProt Entry: [Q6GYQ0](#)

Application Information

Molecular Weight: Predicted, 230 kDa; observed, 260 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB335

Species Reactivity: Human

Applications Tested: Western blotting (WB), flow cytometry (FCM)

Immunogen

Recombinant protein of human RALGAPA1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:1,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Note: This product is for research use only.

Validation Data



Copyright ©2025 Genuin Biotechnologies LLC

RALGAPA1-Alexa Fluor® 647

Flow cytometric analysis of RALGAPA1 expression in HeLa cells using anti-RALGAPA1 antibody (Cat#5246, 1:2,000). Green, isotype control; red, RALGAPA1.


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-RALGAPA1 antibody (Cat#5246). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-RALGAPA1 antibody (Cat#5246, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM