

Catalog #: 52502

Aliases

TMEM87A; Transmembrane Protein 87A; ELKIN1; DKFZP564G2022; Elkin1

Background

Gene Name: TMEM87A

NCBI Gene Entry: [25963](#)

UniProt Entry: [Q8NBN3](#)

Application Information

Molecular Weight: Predicted, 63 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human TMEM87A

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH7.4) containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

