

Catalog #: 52510

Aliases

Ly6g; Ly-6G; Ly-6G.1; Lymphocyte antigen 6G

Background

Gene Name: Ly6g

NCBI Gene Entry: [546644](#)

UniProt Entry: [P35461](#)

Application Information

Molecular Weight: Predicted, 14 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB1470

Species Reactivity: Mouse

Applications Tested: Immunohistochemistry (IHC), flow cytometry (FCM)

Immunogen

Recombinant protein of mouse Ly6g

Isotype

Rabbit IgG

Storage Buffer

Supplied in 50mM Tris-Glycine (pH 7.4) containing 40% glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Immunohistochemistry (IHC): 1:50-1:200

Flow Cytometry (FCM): 1:100-1:200

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

