

Catalog #: 52530

Aliases

ENTPD1; Ectonucleoside Triphosphate Diphosphohydrolase 1; ATPDase; NTPDase-1; SPG64; CD39; Nucleoside Triphosphate Diphosphohydrolase 1; Lymphoid Cell Activation Antigen; Ecto-ATP Diphosphohydrolase 1; ATP Diphosphohydrolase; Ecto-ATPDase 1; Ecto-ATPase 1; CD39 Antigen; Ecto-Apyrase; EC 3.6.1.5; NTPDase1; ATP-DPH; EC 3.6.1

Background

Gene Name: ENTPD1

NCBI Gene Entry: [953](#)

UniProt Entry: [P49961](#)

Application Information

Molecular Weight: Predicted, 58 kDa; observed, 56 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ENTPD1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of HeLa cells using ENTPD1 antibody. Secondary antibody was diluted at 1:20,000.