

Catalog #: 52533

Aliases

TBX21; T-Box Transcription Factor 21; TBLYM; T-Bet; T-Cell-Specific T-Box Transcription Factor T-Bet; T-Box Transcription Factor Expressed In T Cells; T-Box Transcription Factor TBX21; Transcription Factor TBLYM; T-Box Protein 21; T-Box 21; TBET; T-Box Expressed In T Cells; IMD88; T-PET

Background

Gene Name: TBX21

NCBI Gene Entry: [30009](#)

UniProt Entry: [Q9UL17](#)

Application Information

Molecular Weight: Predicted, 58 kDa; observed, 62 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 25GB1495

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human TBX21

Isotype

Mouse IgG2a,Kappa

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.05% BSA and 0.05% Proclin 300.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunocytochemistry (IC): 1:50-1:200

Note: This product is for research use only.

Validation Data

