

Catalog #: 52546

Aliases

TNFAIP8L2; TNF Alpha Induced Protein 8 Like 2; Tumor Necrosis Factor Alpha-Induced Protein 8-Like Protein 2; TNF Alpha-Induced Protein 8-Like Protein 2; Inflammation Factor Protein 20; TNFAIP8-Like Protein 2; FLJ23467; TIPE2; Tumor Necrosis Factor, Alpha-Induced Protein 8-Like Protein 2; Tumor Necrosis Factor, Alpha-Induced Protein 8-Like 2; Tumor Necrosis Factor, Alpha Induced Protein 8 Like 2; Inflammation Factor 20

Background

Gene Name: TNFAIP8L2

NCBI Gene Entry: [Q6P589](#)

UniProt Entry: [79626](#)

Application Information

Molecular Weight: Predicted, 20 kDa; observed, 20 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human TIPE2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blot analysis of the lysates from HeLa cells using TNFAIP8L2 antibody diluted at 1:500.