

Catalog #: 52547

Aliases

ACTR1A; Actin Related Protein 1A; ARP1; Arp1A; ARP1 Actin-Related Protein 1 Homolog A, Centractin Alpha; ARP1 Actin Related Protein 1 Homolog A; Centrosome-Associated Actin Homolog; Alpha-Centractin; Centractin; Actin-RPV; CTRN1; ARP1 (Actin-Related Protein 1, Yeast) Homolog A (Centractin Alpha); Epididymis Secretory Sperm Binding Protein; Centractin Alpha

Background

Gene Name: ACTR1A

NCBI Gene Entry: [10121](#)

UniProt Entry: [P61163](#)

Application Information

Molecular Weight: Predicted, 43 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human ACTZ

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

