

Catalog #: 52548

Aliases

MOGS; Mannosyl-Oligosaccharide Glucosidase; GCS1; Processing A-Glucosidase I ; CWH41; DER7; Glucosidase I; EC 3.2.1.106; Epididymis Secretory Sperm Binding Protein; CDG2B

Background

Gene Name: MOGS

NCBI Gene Entry: [7841](#)

UniProt Entry: [Q13724](#)

Application Information

Molecular Weight: Predicted, 92 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human MOGS

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

