

Catalog #: 52587

Aliases

PF4; Platelet Factor 4; SCYB4; CXCL4; Chemokine (C-X-C Motif) Ligand 4; C-X-C Motif Chemokine 4; Oncostatin-A; Iroplact; PF-4

Background

Gene Name: PF4

NCBI Gene Entry: [5196](#)

UniProt Entry: [P02776](#)

Application Information

Molecular Weight: Predicted, 11 kDa; observed, 13 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human PF4

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PF4 expression in HCT116 (A), mouse brain (B) whole cell lysates.