

Catalog #: 52590

Aliases

DDX49; DEAD-Box Helicase 49; Dbp8; DEAD (Asp-Glu-Ala-Asp) Box Polypeptide 49; Probable ATP-Dependent RNA Helicase DDX49; DEAD Box Protein 49; FLJ10432; EC 3.6.4.13; R27090_2; EC 3.6.1

Background

Gene Name: DDX49

NCBI Gene Entry: [54555](#)

UniProt Entry: [Q9Y6V7](#)

Application Information

Molecular Weight: Predicted, 54 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human DDX49

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

