

**Catalog #: 52599**

## Aliases

IL17RA; Interleukin 17 Receptor A; IL-17RA; CDw217; HIL-17R; CD217; IL17R;  
Interleukin-17 Receptor A; IL-17 Receptor A; Interleukin 17 Receptor; CD217 Antigen;  
CANDF5; IMD51

## Background

Gene Name: IL17RA

NCBI Gene Entry: [23765](#)

UniProt Entry: [Q96F46](#)

## Application Information

Molecular Weight: Predicted, 90 kDa, observed, 90 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

## Immunogen

A synthesized peptide derived from human IL17RA

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

**Note:** This product is for research use only.

## Validation Data

---


### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis of the lysates from HeLa cells using IL17RA antibody.