

Catalog #: 52638

Aliases

DLST; Dihydrolipoamide S-Succinyltransferase; OGDC-E2; KGD2; DLTS; Dihydrolipoyllysine-Residue Succinyltransferase Component Of 2-Oxoglutarate Dehydrogenase Complex, Mitochondrial; Dihydrolipoamide S-Succinyltransferase (E2 Component Of 2-Oxo-Glutarate Complex); 2-Oxoglutarate Dehydrogenase Complex Component E2; EC 2.3.1.61; E2K; Dihydrolipoamide Succinyltransferase Component Of 2-Oxoglutarate Dehydrogenase Complex; Dihydrolipoyllysine-Residue Succinyltransferase; EC 2.3.1; PPGL7; PGL7

Background

Gene Name: DLST

NCBI Gene Entry: [1743](#)

UniProt Entry: [P36957](#)

Application Information

Molecular Weight: Predicted, 49 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human ODO2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry (IHC): 1:50-1:200

Note: This product is for research use only.

Validation Data

