

Catalog #: 52642

Aliases

WNT2B; Wnt Family Member 2B; WNT13; Wingless-Type MMTV Integration Site Family, Member 13; Wingless-Type MMTV Integration Site Family, Member 2B; XWNT2, Xenopus, Homolog Of; Protein Wnt-2b; XWNT2; Protein Wnt-13

Background

Gene Name: WNT2B

NCBI Gene Entry: [7482](#)

UniProt Entry: [Q93097](#)

Application Information

Molecular Weight: Predicted, 44 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB), immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human WNT2B

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry(IHC): 1:50-1:200

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

