

Catalog #: 52646

Aliases

ITGA8; Integrin Subunit Alpha 8; Integrin, Alpha 8; Integrin Alpha-8

Background

Gene Name: ITGA8

NCBI Gene Entry: [8516](#)

UniProt Entry: [P53708](#)

Application Information

Molecular Weight: Predicted, 117 kDa; observed, 130 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB), immunohistochemistry-Paraffin (IHC-P)

Immunogen

A synthesized peptide derived from human ITGA8

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, and 0.02% Proclin 300.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:200-1:1,000

Immunohistochemistry-Paraffin (IHC-P): 1:100-1:500

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Western blotting analysis using anti-ITGA8 antibody (Cat#52646, 1:2,000). Total cell lysates (20 µg) from various cell lines were loaded and separated by SDS-PAGE.