

## Catalog #: 52710

### Aliases

HS1BP3; HCLS1 Binding Protein 3; HCLS1-Binding Protein 3; HS1-Binding Protein 3; HSP1BP-3; HS1-BP3,FLJ14249; HS1-BP3; ETM2

### Background

Gene Name: HS1BP3

NCBI Gene Entry: [64342](#)

UniProt Entry: [Q53T59](#)

### Application Information

Molecular Weight: Predicted, 43 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

### Immunogen

A synthesized peptide derived from human HS1BP3

### Isotype

Rabbit IgG

### Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

### Storage

Store at -20 °C for one year.

### Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

**Note:** This product is for research use only.

### Validation Data

---

#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

