

Catalog #: 52715

Aliases

ANGPTL2; Angiopoietin Like 2; ARP2; HARP; Angiopoietin-Related Protein 2; Angiopoietin-Like Protein 2; Angiopoietin-Like 2

Background

Gene Name: ANGPTL2

NCBI Gene Entry: [23452](#)

UniProt Entry: [Q9UKU9](#)

Application Information

Molecular Weight: Predicted, 57 kDa; observed, 65 kDa

Clonality: Mouse oligoclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human ANGPTL2

Isotype

Mouse IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.09% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-ANGPTL2 antibody (Cat#52715). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-ANGPTL2 antibody (Cat#52715, 1:2,500) and HRP-conjugated goat anti-mouse secondary

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

antibody (Cat#101, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM