

Catalog #: 52728

Aliases

PPIH; Peptidylprolyl Isomerase H; USA-CYP; CYPH; Small Nuclear Ribonucleoprotein Particle-Specific Cyclophilin H; Peptidyl-Prolyl Cis-Trans Isomerase H; Rotamase H; SnuCyp-20; PPIase H; CYP-20; U-SnRNP-Associated Cyclophilin SnuCyp-20; U-SnRNP-Associated Cyclophilin SnuCyp-20; USA-CyP SnuCyp-20; Cyclophilin H; EC 5.2.1.8; MGC5016; Peptidyl Prolyl Isomerase H (Cyclophilin H); CYP20; CypH

Background

Gene Name: PPIH

NCBI Gene Entry: [10465](#)

UniProt Entry: [O43447](#)

Application Information

Molecular Weight: Predicted, 19 kDa; observed, 19 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human PPIH

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-PPIH antibody (Cat#52728). Total lysates (30 µg) were loaded and separated by SDS-PAGE. The blot was incubated with anti-PPIH antibody

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

(Cat#52728, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).