

Catalog #: 52731

Aliases

GAL3ST1; Galactose-3-O-Sulfotransferase 1; CST; Galactosylceramide Sulfotransferase; 3'-Phosphoadenosine-5'-Phosphosulfate:GalCer Sulfotransferase; Cerebroside (3'-Phosphoadenylylsulfate:Galactosylceramide 3') Sulfotransferase; Cerebroside (3'-Phosphoadenylylsulfate:Galactosylceramide 3') Sulfotransferase; 3'-Phosphoadenylylsulfate:Galactosylceramide 3'-Sulfotransferase; Cerebroside Sulfotransferase; GalCer Sulfotransferase; EC 2.8.2.11

Background

Gene Name: GAL3ST1

NCBI Gene Entry: [9514](#)

UniProt Entry: [Q99999](#)

Application Information

Molecular Weight: Predicted, 49 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human GAL3ST1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,000

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Validation Data

