

Catalog #: 52740

Aliases

SLC25A18; Solute Carrier Family 25 Member 18; Solute Carrier Family 25 (Mitochondrial Carrier), Member 18; Solute Carrier Family 25 (Glutamate Carrier), Member 18; Mitochondrial Glutamate Carrier; Glutamate/H(+) Symporter; GC2; GC-2

Background

Gene Name: SLC25A18

NCBI Gene Entry: [83733](#)

UniProt Entry: [Q9H1K4](#)

Application Information

Molecular Weight: Predicted, 34 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human SLC25A18

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

