

Catalog #: 52759

Aliases

CPT1A; Carnitine Palmitoyltransferase 1A; CPT1-L; L-CPT1; CPT1; Carnitine O-Palmitoyltransferase 1, Liver Isoform; Carnitine Palmitoyltransferase 1A (Liver); Succinyltransferase CPT1A; EC 2.3.1.21; CPTI-L; CPT I; Carnitine O-Palmitoyltransferase I, Liver Isoform; Carnitine Palmitoyltransferase I, Liver; EC 2.3.1.-; EC 2.3.1

Background

Gene Name: CPT1A

NCBI Gene Entry: [1374](#)

UniProt Entry: [P50416](#)

Application Information

Molecular Weight: Predicted, 88 kDa; observed, 85 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Mouse, Rat

Applications Tested: Western Blotting (WB)

Immunogen

A synthesized peptide derived from human CPT1A

Isotype

Rabbit IgG

Storage Buffer

Supplied in 0.01M TBS (pH7.4) containing 50% glycerol, 1% BSA and 0.02% Proclin300.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

