

Catalog #: 52761

Aliases

CIDEc; Cell Death Inducing DFFA Like Effector C; CIDE-3; Fsp27; Cell Death Activator CIDE-3; Lipid Transferase CIDEc; FLJ20871; Cell Death-Inducing DFFA-Like Effector Protein C; Cell Death-Inducing DFFA-Like Effector C; Fat-Specific Protein FSP27 Homolog; Fat Specific Protein 27; CIDE3; FPLD5; FSP27

Background

Gene Name: CIDEc

NCBI Gene Entry: [63924](#)

UniProt Entry: [Q96AQ7](#)

Application Information

Molecular Weight: Predicted, 27 kDa; observed, 26 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, Mouse, Rat

Applications Tested: Western Blotting (WB), Immunohistochemistry (IHC), Immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human CIDEc

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 50% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Immunohistochemistry (IHC): 1:50-1:200

Immunocytochemistry (IC): 1:50-1:200

Note: This product is for research use only.

Validation Data

