

Catalog #: 52765

Aliases

HNF4G; Hepatocyte Nuclear Factor 4 Gamma; NR2A2; Nuclear Receptor Subfamily 2 Group A Member; Hepatocyte Nuclear Factor 4-Gamma; HNF-4-Gamma; NR2A3

Background

Gene Name: HNF4G

NCBI Gene Entry: [3174](#)

UniProt Entry: [Q14541](#)

Application Information

Molecular Weight: Predicted, 46 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human HNF4G

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

