

Catalog #: 52766

Aliases

MGAT3; Beta-1,4-Mannosyl-Glycoprotein 4-Beta-N-Acetylglucosaminyltransferase; GNT-III; N-Glycosyl-Oligosaccharide-Glycoprotein N-Acetylglucosaminyltransferase III; Mannosyl (Beta-1,4-)-Glycoprotein Beta-1,4-N-Acetylglucosaminyltransferase; N-Acetylglucosaminyltransferase III; GlcNAc-T III; EC 2.4.1.144; GGNT3; GNT3

Background

Gene Name: MGAT3

NCBI Gene Entry: [4248](#)

UniProt Entry: [Q09327](#)

Application Information

Molecular Weight: Predicted, 61 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human MGAT3

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

