

Anti-Tyrosine Hydroxylase Rabbit Polyclonal Antibody



Catalog #: 52775

Aliases

TH; Tyrosine Hydroxylase; Tyrosine 3-Monooxygenase; DYT5b; Tyrosine 3-Hydroxylase; EC 1.14.16.2; TYH; Dystonia 14; EC 1.14.16; DYT14

Background

Gene Name: TH

NCBI Gene Entry: [7054](#)

UniProt Entry: [P07101](#)

Application Information

Molecular Weight: Predicted, 59 kDa; observed, 60 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB), immunohistochemistry (IHC), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Tyrosine Hydroxylase

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry (IHC): 1:100-1:300

Immunocytochemistry (IC): 1:50-1:200

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

