

Catalog #: 52779

Aliases

NMUR2; Neuromedin U Receptor 2; NMU2R; G-Protein Coupled Receptor TGR-1; G-Protein Coupled Receptor FM-4; Neuromedin-U Receptor 2; NMU-R2; TGR1; Growth Hormone Secretagogue Receptor Family, Member 4; TGR-1; FM-4; FM4

Background

Gene Name: NMUR2

NCBI Gene Entry: [56923](#)

UniProt Entry: [Q9GZQ4](#)

Application Information

Molecular Weight: Predicted, 48 kDa; observed, 46 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human NMUR2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

