

Catalog #: 52784

Aliases

PLAT; Plasminogen Activator, Tissue Type; Tissue-Type Plasminogen Activator; T-Plasminogen Activator; EC 3.4.21.68; Alteplase; Reteplase; T-PA; TPA; Plasminogen Activator, Tissue; Plasminogen/Activator Kringle; EC 3.4.21

Background

Gene Name: PLAT

NCBI Gene Entry: [5327](#)

UniProt Entry: [P00750](#)

Application Information

Molecular Weight: Predicted, 63 kDa; observed, 63 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB), immunohistochemistry (IHC), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human PLAT

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry (IHC): 1:50-1:300

Immunocytochemistry (IC): 1:200-1:1,000

Note: This product is for research use only.

SUPPORT


SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Validation Data



Western blotting analysis of A549 cells using PLAT antibody. Secondary antibody was diluted at 1:20,000.