

Anti-Pan Methylated Lysine Mouse Monoclonal Antibody



Catalog #: 52787

Aliases

Dimethyllysine; Methyllysine; Trimethyllysine; Dimethyl-lysine; Methyl-lysine; Trimethyl-lysine; N epsilon dimethyl lysine; N epsilon mono methyllysine; N epsilon trimethyl lysine; N(epsilon) trimethyl lysine; di methyl-lysine, mono methyl-lysine, Trimethyl-lysine; methyl K pan; di methyl-K, mono methyl-K, Trimethyl-K

Background

NCBI Gene Entry: [56-87-1\(CAS\)](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 25GB6960

Species Reactivity: Species independent

Applications Tested: Western blotting (WB), immunohistochemistry (IHC), immunocytochemistry (IC)

Immunogen

Methylated Lysine

Isotype

Mouse IgG1

Storage Buffer

Supplied in 0.01M TBS (pH7.4) containing 50% glycerol, 1% BSA and 0.02% Proclin300.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Immunohistochemistry (IHC): 1:100-1:500

Immunocytochemistry (IC): 1:100-1:500

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

