

Anti-Cleaved-Caspase-1 p20 Rabbit Polyclonal Antibody



Catalog #: 52790

Aliases

CASP1; Caspase 1; ICE; Caspase-1; IL1BC; Caspase 1, Apoptosis-Related Cysteine Peptidase; Interleukin 1, Beta, Convertase; IL-1 Beta-Converting Enzyme; EC 3.4.22.36; P45; Caspase 1, Apoptosis-Related Cysteine Peptidase (Interleukin 1, Beta, Convertase); Caspase 1, Apoptosis-Related Cysteine Protease (Interleukin 1, Beta, Convertase); Interleukin-1 Beta-Converting Enzyme; Interleukin 1-B Converting Enzyme; Interleukin-1 Beta Convertase; IL1B-Convertase; EC 3.4.22; IL1BCE; CASP-1; IL-1BC

Background

Gene Name: CASP1

NCBI Gene Entry: [834](#)

UniProt Entry: [P29466](#)

Application Information

Molecular Weight: Predicted, 20 kDa; observed, 25 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Caspase-1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Validation Data

