

## Catalog #: 52793

### Aliases

ADRA2A; Adrenoceptor Alpha 2A; ADRAR; ADRA2R; Alpha-2 Adrenergic Receptor Subtype C10; Adrenergic, Alpha-2A-, Receptor; Alpha-2A Adrenergic Receptor; Alpha-2A Adrenoreceptor; Alpha-2AAR Subtype C10; Alpha-2A Adrenoceptor; ADRA2; Alpha-2-Adrenergic Receptor, Platelet Type; Alpha-2A-Adrenergic Receptor; Alpha-2AAR; ALPHA2AAR; FPLD8; ZNF32

### Background

Gene Name: ADRA2A

NCBI Gene Entry: [150](#)

UniProt Entry: [P08913](#)

### Application Information

Molecular Weight: Predicted, 51 kDa; observed, 50 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunocytochemistry (IC)

### Immunogen

A synthesized peptide derived from human ADRA2A

### Isotype

Rabbit IgG

### Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

### Storage

Store at -20 °C for one year.


### Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunocytochemistry (IC): 1:200-1:1,000

**Note:** This product is for research use only.

## Validation Data



Western blotting analysis of JK cells using ADRA2A antibody.