

Catalog #: 52794

Aliases

GPX3; Glutathione Peroxidase 3; Extracellular Glutathione Peroxidase; Glutathione Peroxidase 3 (Plasma); Plasma Glutathione Peroxidase; Selenoprotein GPX3; EC 1.11.1.9; GSHPx-3; GSHPx-P; GPx-P; GPx-3; EC 1.11.1; GPXP

Background

Gene Name: GPX3

NCBI Gene Entry: [2878](#)

UniProt Entry: [P22352](#)

Application Information

Molecular Weight: Predicted, 26 kDa; observed, 24 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human GPX3

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of lysates from A431 cells, primary antibody was diluted at 1:1,000, 4°C over night.