

Catalog #: 52796

Aliases

MANEA; Mannosidase Endo-Alpha; Mandaselin; Glycoprotein Endo-Alpha-1,2-Mannosidase; Endomannosidase; Endo-Alpha Mannosidase; EC 3.2.1.130; FLJ12838; HEndo; Alpha 1,2-Endomannosidase; Mannosidase, Endo-Alpha; EC 3.2.1; ENDO

Background

Gene Name: MANEA

NCBI Gene Entry: [79694](#)

UniProt Entry: [Q5SRI9](#)

Application Information

Molecular Weight: Predicted, 54 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, pig

Applications Tested: Immunohistochemistry (IHC), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human MANEA

Isotype

Rabbit IgG

Storage Buffer

Supplied in 0.01M TBS (pH7.4) containing 50% glycerol, 1% BSA and 0.02% Proclin300.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Immunohistochemistry (IHC): 1:100-1:500

Immunocytochemistry (IC): 1:100-1:500

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

