

**Catalog #: 52800**

## Aliases

CUTC; CutC Copper Transporter; CGI-32; Copper Homeostasis Protein CutC Homolog; CutC Copper Transporter Homolog (E. Coli); CutC Copper Transporter Homolog

## Background

Gene Name: CUTC

NCBI Gene Entry: [51076](#)

UniProt Entry: [Q9NTM9](#)

## Application Information

Molecular Weight: Predicted, 29 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, Mouse

Applications Tested: Western Blotting (WB)

## Immunogen

A synthesized peptide derived from human CUTC

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

