

Catalog #: 52801

Aliases

PTH LH; Parathyroid Hormone Like Hormone; PTHRP; PLP; PTHR; HHM; Parathyroid Hormone-Like Hormone Preproprotein; Parathyroid Hormone-Related Protein; Osteostatin; PTH-RP; Parathyroid Hormone-Related Protein Preproprotein; Parathyroid Hormone-Like Related Protein; Parathyroid Hormone-Like Protein; PTH-Related Protein; PTHrP; BDE2

Background

Gene Name: PTH LH

NCBI Gene Entry: [5744](#)

UniProt Entry: [P12272](#)

Application Information

Molecular Weight: Predicted, 20 kDa; observed, 20 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, Mouse, Rat

Applications Tested: Western Blotting (WB), Immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human PTH LH

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry (IHC): 1:50-1:300

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

