

Catalog #: 52803

Aliases

FABP7; Fatty Acid Binding Protein 7; B-FABP; BLBP; Brain-Type Fatty Acid-Binding Protein; Fatty Acid-Binding Protein, Brain; Brain Lipid-Binding Protein; FABPB; MRG; Mammary-Derived Growth Inhibitor-Related; Mammary-Derived Growth Inhibitor Related; Fatty Acid Binding Protein 7, Brain; Hypothetical Protein DKFZp547J2313; Fatty Acid-Binding Protein 7; Brain Lipid Binding Protein

Background

Gene Name: FABP7

NCBI Gene Entry: [2173](#)

UniProt Entry: [O15540](#)

Application Information

Molecular Weight: Predicted, 15 kDa; observed, 13 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB6970

Species Reactivity: Human, Mouse

Applications Tested: Western Blotting (WB), Immunohistochemistry (IHC), Immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human BLBP

Isotype

Rabbit IgG

Storage Buffer

Supplied in 0.01M TBS (pH7.4) containing 50% glycerol, 1% BSA and 0.02% Proclin300.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry (IHC): 1:100-1:500

Immunocytochemistry (IC): 1:50

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Note: This product is for research use only.

Validation Data

