

Catalog #: 52805

Aliases

TNF; Tumor Necrosis Factor; TNF-Alpha; TNFSF2; TNFA; DIF; Tumor Necrosis Factor Ligand Superfamily Member 2; TNF-A 3 4Tumor Necrosis Factor (TNF Superfamily, Member 2); Tumor Necrosis Factor Ligand 1F; Tumor Necrosis Factor-Alpha; Tumor Necrotic Factor Alpha; TNF Superfamily, Member 2; TNF, Macrophage-Derived; TNF, Monocyte-Derived; APC1 Protein; Cachectin; IMD127; TNLG1F

Background

Gene Name: TNF

NCBI Gene Entry: [7124](#)

UniProt Entry: [P01375](#)

Application Information

Molecular Weight: Predicted, 26 kDa; observed, 18 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB6975

Species Reactivity: Human, Mouse, Rat

Applications Tested: Western Blotting (WB), Immunohistochemistry (IHC),

Immunocytochemistry (IC), Immunoprecipitation (IP)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.05% Proclin 300 and 0.05% BSA.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:2,000-1:10,000

Immunohistochemistry (IHC): 1:2,000-1:10,000

Immunocytochemistry (IC): 1:200-1:1,000

Immunoprecipitation (IP): 1:50-1:200

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Validation Data

