

Anti-SARS-COV-2 Spike Protein S2 Rabbit Polyclonal Antibody



Catalog #: 52807

Aliases

SARS-COV-2 Spike Protein S2; Spike glycoprotein; S glycoprotein; E2; Peplomer protein

Background

UniProt Entry: [P0DTC2](#)

Application Information

Molecular Weight: Predicted, 141 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Humanvirus

Applications Tested: Western Blotting (WB)

Immunogen

A synthesized peptide derived from SARS-COV-2 Spike Protein S2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

