

Catalog #: 52808

Aliases

IL17A; Interleukin 17A; IL-17A; IL-17; CTLA8; IL17; Interleukin 17 (Cytotoxic T-Lymphocyte-Associated Serine Esterase 8); Cytotoxic T-Lymphocyte-Associated Protein 8; Cytotoxic T-Lymphocyte-Associated Antigen 8; Interleukin-17A; CTLA-8; Interleukin-17; ILA17

Background

Gene Name: IL17A

NCBI Gene Entry: [3605](#)

UniProt Entry: [Q16552](#)

Application Information

Molecular Weight: Predicted, 18 kDa; observed, 22 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Mouse, Rat

Applications Tested: Western Blotting (WB), Immunohistochemistry (IHC), Immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human IL-17

Isotype

Rabbit IgG

Storage Buffer

Supplied in 0.01M TBS (pH7.4) containing 50% glycerol, 1% BSA and 0.02% Proclin300.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry (IHC): 1:100-1:500

Immunocytochemistry (IC): 1:100-1:500

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Validation Data

