

Catalog #: 52810

Aliases

ITK; IL2 Inducible T Cell Kinase; EMT; LYK; PSCTK2; Interleukin-2-Inducible T-Cell Kinase; Tyrosine-Protein Kinase ITK/TSK; IL-2-Inducible T-Cell Kinase; T-Cell-Specific Kinase; EC 2.7.10.2; Kinase EMT; Homolog Of Mouse T-Cell Itk/Tsk; IL2-Inducible T-Cell Kinase; IL2 Inducible T-Cell Kinase; Tyrosine-Protein Kinase LYK; Tyrosine-Protein Kinase Lyk; EC 2.7.10; LPFS1

Background

Gene Name: ITK

NCBI Gene Entry: [3702](#)

UniProt Entry: [Q08881](#)

Application Information

Molecular Weight: Predicted, 72 kDa; observed, 72 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, Mouse

Applications Tested: Western Blotting (WB), Immunohistochemistry (IHC), Immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human ITK

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry (IHC): 1:50-1:200

Immunocytochemistry (IC): 1:50-1:200

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Note: This product is for research use only.

Validation Data

