

Catalog #: 52811

Aliases

SPHK1; Sphingosine Kinase 1; SPHK; Acetyltransferase SPHK1; SPK 1; SK 1; EC 2.7.1.91; EC 2.3.1.-; SK1; SPK

Background

Gene Name: SPHK1

NCBI Gene Entry: [8877](#)

UniProt Entry: [Q9NYA1](#)

Application Information

Molecular Weight: Predicted, 43 kDa; observed, 42 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, Mouse, Rat

Applications Tested: Western Blotting (WB)

Immunogen

A synthesized peptide derived from human SPHK1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

