

Catalog #: 52815

Aliases

FANCM; FA Complementation Group M; KIAA1596; FAAP250; Fanconi Anemia-Associated Polypeptide Of 250 KDa; Fanconi Anemia Complementation Group M; ATP-Dependent RNA Helicase FANCM; Fanconi Anemia Group M Protein; Protein Hef Ortholog; Protein FACM; EC 3.6.4.13; EC 3.6.1; SPGF28; POF15

Background

Gene Name: FANCM

NCBI Gene Entry: [57697](#)

UniProt Entry: [Q8IYD8](#)

Application Information

Molecular Weight: Predicted, 232 kDa; observed, 225 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, Mouse

Applications Tested: Western Blotting (WB)

Immunogen

A synthesized peptide derived from human FANCM

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

