

Anti-Transferrin Receptor Rabbit Monoclonal Antibody



Catalog #: 52825

Aliases

TFRC; Transferrin Receptor; P90; CD71; TFR1; Transferrin Receptor Protein 1; T9; TR; Transferrin Receptor (P90, CD71); CD71 Antigen; IMD46; TRFR; TfR1; Trfr; TFR; TfR

Background

Gene Name: TFRC

NCBI Gene Entry: [7037](#)

UniProt Entry: [P02786](#)

Application Information

Molecular Weight: Predicted, 85 kDa; observed, 84 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB7115

Species Reactivity: Human, Mouse, Rat

Applications Tested: Western Blotting (WB), Immunohistochemistry (IHC), Immunocytochemistry (IC), Immunoprecipitation (IP)

Isotype

Rabbit IgG,Kappa

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.05% Proclin 300.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Immunohistochemistry (IHC): 1:200-1:1,000

Immunocytochemistry (IC): 1:50-1:200

Immunoprecipitation (IP): 1:50-1:200

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

