

Catalog #: 52838

Aliases

RPS27A; Ribosomal Protein S27a; Ubiquitin Carboxyl Extension Protein 80; UBCEP80; Uba80; S27A; ES31; Ubiquitin-Ribosomal Protein ES31 Fusion Protein; UBCEP1; Ubiquitin And Ribosomal Protein S27a; Ubiquitin-40S Ribosomal Protein S27a; Epididymis Luminal Protein 112; 40S Ribosomal Protein S27a; Ubiquitin-CEP80; Ubiquitin C; HEL112; CEP80; UBA80; UBC

Background

Gene Name: RPS27A

NCBI Gene Entry: [6233](#)

UniProt Entry: [P62979](#)

Application Information

Molecular Weight: Predicted, 18 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human RPS27A

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

