

Catalog #: 5284

Aliases

TECPR1; Tectonin Beta-Propeller Repeat Containing 1; KIAA1358; Tectonin Beta-Propeller Repeat-Containing Protein 1; DKFZP434B0335; FLJ23419; FLJ90593

Background

Gene Name: TECPR1

NCBI Gene Entry: [25851](#)

UniProt Entry: [Q7Z6L1](#)

Application Information

Molecular Weight: Predicted, 130 kDa; observed, 135 kDa

Clonality: Mouse oligoclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB), flow cytometry (FCM)

Immunogen

Recombinant protein of human TECPR1

Isotype

Mouse IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:400-1:2,000

Flow Cytometry (FCM): 1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Flow cytometric analysis of TECPR1 expression in HT-1080 cells using anti-TECPR1 antibody (Cat#5284, 1:2,000). Green, isotype control; red, TECPR1.



Western blotting analysis using anti-TECPR1 antibody (Cat#5284). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-TECPR1 antibody (Cat#5284, 1:2,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQTM ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM