

Catalog #: 52840

Aliases

TMOD3; Tropomodulin 3; Ubiquitous Tropomodulin; UTMOD; Tropomodulin-3; Tropomodulin 3 (Ubiquitous); U-Tmod

Background

Gene Name: TMOD3

NCBI Gene Entry: [29766](#)

UniProt Entry: [Q9NYL9](#)

Application Information

Molecular Weight: Predicted, 40 kDa; observed, 38 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human TMOD3

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

