

Catalog #: 52846

Aliases

PLD2; Phospholipase D2; Choline Phosphatase 2; Phosphatidylcholine-Hydrolyzing
Phospholipase D2; EC 3.1.4.4; PLD1C; PLD; HPLD2

Background

Gene Name: PLD2

NCBI Gene Entry: [5338](#)

UniProt Entry: [O14939](#)

Application Information

Molecular Weight: Predicted, 106 kDa; observed, 50 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human PLD2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

Immunohistochemistry (IHC): 1:50-1:300

Note: This product is for research use only.