

Catalog #: 52847

Aliases

IMP4; IMP U3 Small Nucleolar Ribonucleoprotein 4; BXDC4; U3 Small Nucleolar Ribonucleoprotein Protein IMP4; Brix Domain-Containing Protein 4; U3 SnoRNP Protein IMP4; MGC19606; IMP4, U3 Small Nucleolar Ribonucleoprotein, Homolog (Yeast); IMP4, U3 Small Nucleolar Ribonucleoprotein, Homolog; IMP4 Homolog, U3 Small Nucleolar Ribonucleoprotein; U3 SnoRNP Protein 4 Homolog

Background

Gene Name: IMP4

NCBI Gene Entry: [92856](#)

UniProt Entry: [Q96G21](#)

Application Information

Molecular Weight: Predicted, 34 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human IMP4

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

